

Jelly Bean

Long considered the standard in the industry, our famous Jelly Bean switch features a 2.5-in/6.4-cm activation surface and is recommended for users who are able to access a smaller target area.

Switch tops can be removed and replaced with the color of your choice: Red, Blue, Yellow, or Green. Also includes a clear snap cap for symbol use. It's all up to your own style, tastes, or specific vision needs.

FOR MORE INFORMATION

If you think this device may be beneficial to you, please contact your SSA and request an individual consult. If you are a provider seeking training or consultation contact us at (419) 782-6621.



PRODUCT FEATURES

Wired Connection

Pressure Activation

Auditory Click, Tactile Feedback

All Ages

