



Defiance DD

County Board of Developmental Disabilities

UNITE. INSPIRE. BECOME.

Bintiva Peanut Ball

Our peanut ball is crafted from professional grade materials, and is similar to a standard stability ball. The difference being, that, peanut balls only move back and forth in one direction, versus the stability ball that is multidirectional. This provides more stability and control, and makes it ideal for therapy and rehabilitation use with children or seniors. The unique shape is also favored by midwives, and doulas, for use as a birthing ball.

FOR MORE INFORMATION

If you think this device may be beneficial to you, please contact your SSA and request an individual consult. If you are a provider seeking training or consultation contact us at (419) 782-6621.



PRODUCT FEATURES

Polyvinyl Chloride

3lbs.

Durable and Longlasting

Anti Burst

Versatile

